

 KEWAUNEE®

ADJUSTAbench™



THE POWER OF INGENUITY



# ADJUSTAbench™

ADJUSTAbench™ answers the need for optimized ergonomics in the assembly plant, the laboratory, the clean room —

wherever productivity goes up and absence rates go down in direct proportion to worker comfort.


As the name implies, ADJUSTAbench permits adjusting the worksurface height — at the push of a button — to any dimension

within its **28" to 44" design range**. Smoothly and quietly assuming the technician's precise comfort setting, this unique furniture solution makes particular sense in multi-task or multi-user work environments.

Select from standard widths of **36", 48", 60", or 72"**, or **specify a custom workstation for your unique needs**. Choose from an extensive list of accessories to configure the unit for its particular assignment — shelving, powered back-splashes, utility raceways, swing arms, task lighting, hanging drawers, foot rests, tool bars, and more. And you can specify an ADJUSTAbench in a wide array of colors and laminate options.

If static's a concern, **specify the ESD worktop for total electrostatic discharge protection**. It comes complete with two grounding terminuses — one serving as the station's single grounding point, the other for continuity testing.

ADJUSTAbench is another example of **The Power of Ingenuity** found at Kewaunee, providing an efficient solution for maximizing productivity and worker comfort.



*Shelf and worktop move as one, maintaining an ergonomic relationship. Intelligent column design provides infinite adjustment of accessories above the worksurface. Rugged, rectangular steel framing under the high-density laminate worksurface supports up to 1,000 pounds.*



*ADJUSTAbench can be programmed to recall three specific settings for quick, pre-set repositioning from shift to shift.*



*ADJUSTAbench's telescopic elevation mechanism is shielded within a protective housing. All metal components feature a durable, chemical-resistant finish.*

 **KEWAUNEE®**  
Technical Products Group

THE POWER OF INGENUITY

2700 West Front Street  
Statesville NC 28677  
1-800-824-6626  
Email: [custsvc@kewaunee-tpg.com](mailto:custsvc@kewaunee-tpg.com)