


Hood Enclosures & Finished Backs

Fume Hood Ceiling Enclosures




Fume Hood Ceiling Enclosures are designed to fill the space between the top of the hood and the ceiling to provide a finished appearance. They are available in height from 9" to 30" in 1" increments and in lengths to match

the hoods. **CF** Style enclosures are 1 1/8" deeper than **HV** Style enclosures to accommodate thicker sash assemblies.

All are provided with a removable front panel for easy access to the light fixture.

Enclosure Type	Material	Enclosure Height	Enclosure Depth	Overall Length	Style	Option
CE	M					
CE = Ceiling Enclosure	M = Powder Coated Steel					FB = Finished Back
Enclosure Height						Styles
18 = 18"	22 = 22"	26 = 26"				Horizontal Sash Hoods = HV V07-V67
19 = 19"	23 = 23"	27 = 27"				& Vertical Sash Bench Hoods V05-V10-V15
20 = 20"	24 = 24"	28 = 28"				V25-V30-V40-V45
21 = 21"	25 = 25"	29 = 29"				Combo Sash Hoods = CF V06-V11-V16-V26-V36
		30 = 30"				& Floor Mounted Hoods V65-V66-V90
Enclosure Depth						Hood Overall Length
24 = 24"	30 = 30"	36 = 36"				4'-0" / 48" = 48
						5'-0" / 60" = 60
						6'-0" / 72" = 72
						8'-0" / 96" = 96
						10'-0" / 120" = 20
						12'-0" / 144" = 44

Fume Hood Finished Backs



Fume Hood Finished Backs are furnished in two pieces

Fume Hood Finished Backs are designed to enclose the back of the fume hood when it is exposed to view and are removable without the use of tools for when access is required.

General Purpose Hood 57" High	ADA Hood 60" High	Length
VFBM570148	VFBM600148	4'-0" / 48"
VFBM570160	VFBM600160	5'-0" / 60"
VFBM570172	VFBM600172	6'-0" / 72"
VFBM570196	VFBM600196	8'-0" / 96"
VFBM570120		
VFBM570144		
LX Series Hood 69" High	Floor Mounted Hood 92" High	Length
VFBM690148	VFBM920148	4'-0" / 48"
VFBM690160	VFBM920160	5'-0" / 60"
VFBM690172	VFBM920172	6'-0" / 72"
VFBM690196	VFBM920196	8'-0" / 96"
VFBM690120	VFBM920120	10'-0" / 120"
VFBM690144	VFBM920144	12'-0" / 144"